

H1N1 HA recombinant protein (C-His)



Catalog Number: 514924

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Bel/1942 (H1N1)) hemagglutinin (ABD62843.1) (Met1-Gln529), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H1N1

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 18

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Bel/1942 (H1N1)) consists of 523 amino acids and predicts a molecular mass of 59.3 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % gly, pH 8.0.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.


Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



H1N1 HA Protein 15521



Influenza A H1N1 (A/Bel/1942) Hemagglutinin / HA Protein (His Tag)