

# Itolizumab reference antibody



Catalog Number: 176924

## Product name

Itolizumab reference antibody

## Specificity

Human CD6

## Antibody description

Anti-CD6 Antibody (Itolizumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations



## SDS-PAGE

Anti-CD6 Antibody (itolizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-CD6 Antibody (itolizumab) is 100.00%, determined by SEC-HPLC.