

# Anti-phospho-PML Protein (Ser518) Rabbit antibody



Catalog Number: 167924

## Product name

Anti-phospho-PML Protein (Ser518) Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

phospho-PML Protein (Ser518) Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser518 of human PML Protein

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/1000

IHC: 1/50

## Validations



Western blot detection of phospho-PML Protein (Ser518) in A549,HL-60,U2OS cell lysates using phospho-PML Protein (Ser518) Rabbit mAb(1:1000 diluted).Predicted band size:98kDa.Observed band size:98-117kDa.