

# AK3 antibody

Catalog Number: 107924

## Product name

AK3 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:200-1:2000

IHC: 1:20-1:200

## Antibody description

AK3 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in human placenta tissue, BxPC-3 cells, HEK-293 cells, HepG2 cells, human kidney tissue, L02 cells, MCF7 cells, mouse heart tissue, mouse kidney tissue, Y79 cells. Positive IHC detected in human pancreas cancer tissue, human heart tissue. Observed molecular weight by Western-blot: 26 kDa

## Validations




## Preparation

This antibody was obtained by immunization of AK3 recombinant protein (Accession Number: NM\_001199854). Purification method: Antigen affinity purified.

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107924(AK3 antibody) at dilution of 1:500

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:107924(AK3 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC, IP

## Dilutions

# AK3 antibody

Catalog Number: 107924



IP Result of anti-AK3 (IP:Catalog No:107924, 3ug; Detection:Catalog No:107924 1:500) with mouse kidney tissue lysate 4400ug.