

# CAMKK2 antibody



Catalog Number: 108924

## Product name

CAMKK2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CAMKK2 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive WB detected in PC-3 cells, fetal human brain tissue, HeLa cells, human brain tissue, human heart tissue, Jurkat cells. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of CAMKK2 recombinant protein (Accession Number: NM\_172226). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108924(CAMKK2 antibody) at dilution of 1:500



Immunofluorescent analysis of PC-3 cells using Catalog No:108924(CAMKK2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)