

# DHRS3 antibody

Catalog Number: 109924

## Product name

DHRS3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DHRS3 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Positive IP detected in A375 cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of DHRS3 recombinant protein (Accession Number: NM\_004753). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:109924(DHRS3 antibody) at dilution of 1:500



IP Result of anti-DHRS3 (IP:Catalog No:109924, 4ug; Detection:Catalog No:109924 1:500) with A375 cells lysate 3600ug.