

ALOX12 Polyclonal Antibody



Catalog Number: 162224

Product name

ALOX12 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ALOX12

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 483-663 of human ALOX12 (NP_000688.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - ALOX12 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using ALOX12 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ALOX12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using ALOX12 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ALOX12 Polyclonal Antibody

ALOX12 Polyclonal Antibody



Catalog Number: 162224

Immunohistochemistry of paraffin-embedded

mouse brain using ALOX12 antibody at dilution of 1:100 (40x lens).