

NPTX1 antibody



Catalog Number: 113224

Product name

NPTX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NPTX1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002522). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113224(NPTX1 antibody) at dilution of 1:400



Immunofluorescent analysis of Hela cells, using NPTX1 antibody Catalog No:113224 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).