

# FRS2 antibody

Catalog Number: 107224

## Product name

FRS2 antibody

## Specificity

Human; other species not tested.

## Antibody description

FRS2 Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, HeLa cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of FRS2 recombinant protein (Accession Number: NM\_006654). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG3

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with (FRS2 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:107224(FRS2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)