

Anti-Cystatin C (5C3) Mouse antibody



Catalog Number: 168224

Product name

Anti-Cystatin C (5C3) Mouse antibody

Specificity

Human

Antibody description

Cystatin C (5C3) Mouse monoclonal antibody

Preparation

Antigen: Recombinant human cystatin c protein.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications

ELISA

Validations



Standard Curve for Human Cystatin C: Capture Antibody Mouse mAb [6F12-C7-D8] to Human Cystatin C at 4u03bcg/ml and Detector Antibody Mouse mAb [168224:5C3-C7.200912-3B11:3B11-E9] to Human Cystatin C at 0.5u03bcg/ml.