

# fetuin-B antibody



Catalog Number: 110624

## Product name

fetuin-B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

fetuin-B Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, mouse liver tissue. Positive WB detected in human blood tissue, human liver tissue. Observed molecular weight by Western-blot: 50-60kd

## Preparation

This antibody was obtained by immunization of fetuin-B recombinant protein (Accession Number: NM\_001375587). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:110624(fetuin-B antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human liver using Catalog No:110624(fetuin-B antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110624(fetuin-B antibody) at dilution of 1:50 (under 40x lens)