

MIOS antibody

Catalog Number: 112624

Product name

MIOS antibody

Specificity

Human; other species not tested.

Antibody description

MIOS Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of MIOS recombinant protein (Accession Number: BC005883). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112624(MIOS Antibody) at dilution of 1:600



IP Result of anti-MIOS (IP:Catalog No:112624, 4ug; Detection:Catalog No:112624 1:300) with HeLa cells lysate 4000ug.