

# Anti-PTPRA antibody



Catalog Number: 103624

## Product name

Anti-PTPRA antibody

## Immunogen

[Human PTPRA \(aa 174-793, His & GST Tag\) recombinant protein](#)

## Specificity

Human PTPRA

## Antibody description

Rabbit polyclonal to PTPRA

## Preparation

Produced in rabbits immunized with purified, recombinant Human PTPRA (rh PTPRA; P18433-2; Ala174-Lys793). PTPRA specific IgG was purified by Human PTPRA affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-2 µg/ml

## Validations



PTPRA Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human PTPRA in human hepatoma with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).