

Human ATG10/APG10L/PP12616 (His tag) recombinant protein



Catalog Number: 519624

General Information

Gene Name Synonym

Ubiquitin-Like-Conjugating Enzyme ATG10, Autophagy-Related Protein 10, APG10-Like, ATG10

Protein Construction

A DNA sequence encoding the Human ATG10/APG10L/PP12616 Met1-Thr190 is expressed with a T7 tag at the N-terminus, 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ATG10/APG10L/PP12616 (His tag)
recombinant protein