

ENPP4 antibody



Catalog Number: 110324

Product name

ENPP4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ENPP4 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 52-60kd

Preparation

This antibody was obtained by immunization of ENPP4 recombinant protein (Accession Number: NM_014936). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:110324(ENPP4 antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human brain using Catalog No:110324(ENPP4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:110324(ENPP4 antibody) at dilution of 1:100 (under 40x lens)