

ODF1 antibody



Catalog Number: 113324

Product name

ODF1 antibody

Specificity

Human, Rat, Mouse; other species not tested.

Antibody description

ODF1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in rat testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 28,56 kDa

Preparation

This antibody was obtained by immunization of ODF1 recombinant protein (Accession Number: NM_024410). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113324(ODF1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113324(ODF1 Antibody) at dilution of 1:50 (under 40x lens)