

# PTPN9 antibody



Catalog Number: 114324

## Product name

PTPN9 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PTPN9 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A549 cells, mouse pancreas tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of PTPN9 recombinant protein (Accession Number: NM\_002833). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114324(PTPN9 Antibody) at dilution of 1:1000