

Anti-DENND2B antibody



Catalog Number: 176324

Product name

Anti-DENND2B antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to DENND2B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human st5. aa range:411-460

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations (kD)




Fig1.; Western Blot analysis of various cells using ST5 Polyclonal Antibody