

ATPIF1 antibody

Catalog Number: 108324

Product name

ATPIF1 antibody

IP: 1:200-1:2000

IF: 1:50-1:500

Specificity

Human, Mouse; other species not tested.

Antibody description

ATPIF1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 12 kDa

Validations




Preparation

This antibody was obtained by immunization of ATPIF1 recombinant protein (Accession Number: NM_016311). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108324(ATPIF1 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

IP Result of anti-ATPIF1 (IP:Catalog No:108324, 4ug; Detection:Catalog No:108324 1:500) with HeLa cells lysate 3440ug.

Ig Type

Rabbit IgG


Applications

ELISA, WB, IP, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog

ATPIF1 antibody



Catalog Number: 108324

No:108324(ATPIF1 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)