

Human BCL2L11/BIM (His tag) recombinant protein



Catalog Number: 518324

General Information

Gene Name Synonym

Bcl-2-like protein 11; Bcl2-L-11; Bcl2-interacting mediator of cell death; BCL2L11; BIM; BIML

Protein Construction

A DNA sequence encoding the Human BCL2L11/BIM Met1-Arg120 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM HEPES, 150mM NaCl, 10% glycerol, 2mM DTT, pH8.0

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human BCL2L11/BIM (His tag) recombinant protein