

STK38L antibody



Catalog Number: 115724

Product name

STK38L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

STK38L Rabbit Polyclonal antibody. Positive WB detected in A375 cells, HepG2 cells. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of STK38L recombinant protein (Accession Number: XM_006719058). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115724(STK38L antibody) at dilution of 1:400