

# Anti-TIMP1 antibody

Catalog Number: 176724

## Product name

Anti-TIMP1 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to TIMP1

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*TBS (pH7.4), 1%BSA, 40%Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IP, IF

## Dilutions

WB: 1:1,000

IP: 1-2 µg per 100-500 µg of total protein

IF: 1:50-500

## Validations




Fig1: Western blot analysis of TIMP-1 expression in non-transfected(A) and human TIMP-1 transfected (B) 293T whole cell lysates and human recombinant TIMP-1 (C).