

SCTR Polyclonal Antibody



Catalog Number: 164034

Product name

SCTR Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to SCTR

Preparation

Antigen: Recombinant protein of human SCTR

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Western blot - SCTR Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using SCTR antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.