

# Anti-YTHDC2 Rabbit antibody



Catalog Number: 168034

## Product name

Anti-YTHDC2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

YTHDC2 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human YTHDC2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000

IHC: 1/50

FC: 1/50

IP: 1/20

## Validations



Western blot detection of YTHDC2 in HeLa, A549, HL-60, U251, U87-MG cell lysates using YTHDC2 Rabbit mAb(1:1000 diluted). Predicted band size: 160kDa. Observed band size: 160kDa.