

# TSHZ3 antibody



Catalog Number: 116434

## Product name

TSHZ3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TSHZ3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse uterus tissue. Observed molecular weight by Western-blot: 105-120 kDa

## Preparation

This antibody was obtained by immunization of TSHZ3 recombinant protein (Accession Number: BC127095). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116434(TSHZ3 Antibody) at dilution of 1:600