

CNOT3 antibody

Catalog Number: 109434

Product name

CNOT3 antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CNOT3 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HEK-293 cells, HT-1080 cells, human brain tissue, K-562 cells, mouse brain tissue. Positive IP detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of CNOT3 recombinant protein (Accession Number: NM_014516). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109434(CNOT3 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using CNOT3 antibody Catalog No:109434 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-CNOT3 (IP:Catalog No:109434,

CNOT3 antibody



Catalog Number: 109434

4ug; Detection:Catalog No:109434 1:1000) with HeLa cells lysate 3200ug.