

GALC antibody

Catalog Number: 110834

Product name

GALC antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GALC Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of GALC recombinant protein (Accession Number: BC036518). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110834(GALC antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110834(GALC antibody) at dilution of 1:50 (under 10x lens)