

Anti-TXN antibody



Catalog Number: 142834

Product name

Anti-TXN antibody

Immunogen

[Mouse Thioredoxin/TXN/SASP recombinant protein](#)

Specificity

Mouse Thioredoxin/TXN/SASP

Antibody description

Rabbit Polyclonal to Mouse TXN

Preparation

Produced in rabbits immunized with purified, recombinant Mouse Thioredoxin/TXN/SASP. Thioredoxin/TXN/SASP specific IgG was purified by Mouse Thioredoxin/TXN/SASP affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, ICC/IF, IF

Dilutions


WB: 1-5 µg/ml

ELISA: 0.1-0.2 µg/mL


This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Mouse Thioredoxin/TXN/SASP.

ICC/IF: 0.5-3µg/mL

Validations



Mouse Thioredoxin/TXN/SASP Immunofluorescence(IF) 15679



Mouse Thioredoxin/TXN/SASP Western blot (WB) 15678