

General Information

Gene Name Synonym

Interferon alpha-5; MCG123694

Protein Construction

A DNA sequence encoding the mouse IFNA5 (NP_034635.2) (Cys24-Glu189) was expressed with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Yeast

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Cys 24

Molecular Mass

The recombinant mouse IFNA5 consists of 176

amino acids and predicts a molecular mass of 20.3 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

A hardcopy of COA with reconstitution instruction is sent along with the products. Please refer to it for detailed information.

SDS-PAGE



Mouse IFNA5/IFNaG Protein 14762