

Anti-Phospho-CDK1/2/3 (Thr14) Rabbit antibody



Catalog Number: 168834

Product name

Rabbit IgG

Anti-Phospho-CDK1/2/3 (Thr14) Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

Phospho-CDK1/2/3 (Thr14) Rabbit monoclonal antibody

IHC: 1/50

ICC/IF: 1/50

Preparation

IP: 1/20

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr14 of human Cdk1

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of Phospho-CDK1/CDK2/CDK3 (Thr14) in HeLa, CHO-K1 cell lysates using Phospho-CDK1/CDK2/CDK3 (Thr14) Rabbit mAb(1:500 diluted). Predicted band size:34kDa. Observed band size:34kDa.

