

GPR75 antibody



Catalog Number: 111134

Product name

GPR75 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GPR75 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of GPR75 recombinant protein (Accession Number: NM_006794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111134(GPR75 Antibody) at dilution of 1:100 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111134(GPR75 Antibody) at dilution of 1:100 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).