

CD41/Integrin alpha 2b antibody



Catalog Number: 107134

Product name

CD41/Integrin alpha 2b antibody

Specificity

Human; other species not tested.

Antibody description

CD41/Integrin alpha 2b Mouse Monoclonal antibody. Positive WB detected in human spleen tissue. Observed molecular weight by Western-blot: 125 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000419). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107134(CD41/Integrin alpha 2b Antibody) at dilution of 1:1600