

# Anti-Syndecan-1 Mouse antibody



Catalog Number: 169134

## Product name

Anti-Syndecan-1 Mouse antibody

## Specificity

Human, Mouse

## Antibody description

Syndecan-1 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human SDC1 (AA: 28-171) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

ICC/IF, FC

## Validations




Western blot analysis using SDC1 mAb against



human SDC1 (AA)

Western blot analysis using SDC1 mAb against  
HEK293 (1) and SDC1 (AA)




Immunofluorescence analysis of HepG2 cells using  
SDC1 mouse mAb (green). Blue


# Anti-Syndecan-1 Mouse antibody



Catalog Number: 169134



Flow cytometric analysis of A431 cells using SDC1 mouse mAb (green) and negative control (red).



Western blot analysis using SDC1 mouse mAb against Hela (1), MCF-7 (2), HepG2 (3), T47D (4), Jurkat (5), NIH/3T3 (6) cell lysate.