

Anti-Syndecan-1 Mouse antibody



Catalog Number: 169134

Product name

Anti-Syndecan-1 Mouse antibody

Specificity

Human, Mouse

Antibody description

Syndecan-1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human SDC1 (AA: 28-171) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ICC/IF, FC

Validations




Western blot analysis using SDC1 mAb against



human SDC1 (AA)


Western blot analysis using SDC1 mAb against
HEK293 (1) and SDC1 (AA




Immunofluorescence analysis of HepG2 cells using
SDC1 mouse mAb (green). Blue

Anti-Syndecan-1 Mouse antibody

Catalog Number: 169134



Flow cytometric analysis of A431 cells using SDC1 mouse mAb (green) and negative control (red).



Western blot analysis using SDC1 mouse mAb against HeLa (1), MCF-7 (2), HepG2 (3), T47D (4), Jurkat (5), NIH/3T3 (6) cell lysate.