

FAU antibody

Catalog Number: 110534

Product name

FAU antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAU Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human colon cancer tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 14kd

Preparation

This antibody was obtained by immunization of FAU recombinant protein (Accession Number: NM_001997). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110534(FAU antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:110534(FAU Antibody) at dilution of 1:50 (under 40x lens)