

RAMP2 antibody

Catalog Number: 114534

Product name

RAMP2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAMP2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, SH-SY5Y cells. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 20kd

Preparation

This antibody was obtained by immunization of RAMP2 recombinant protein (Accession Number: NM_005854). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114534(RAMP2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114534(RAMP2 antibody) at dilution of 1:100 (under 10x lens)