

# SPAG6 antibody

Catalog Number: 115534

## Product name

SPAG6 antibody

## Immunogen

[Human SPAG6 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SPAG6 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, Jurkat cells, PC-3 cells. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 50-55 kDa

## Preparation

This antibody was obtained by immunization of SPAG6 recombinant protein (Accession Number: NM\_172242). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115534(SPAG6 Antibody) at dilution of 1:300



IP Result of anti-SPAG6 (IP:Catalog No:115534, 3ug; Detection:Catalog No:115534 1:300) with mouse testis tissue lysate 4000ug.