

UBE2R2 antibody



Catalog Number: 116534

Product name

UBE2R2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBE2R2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, mouse spleen tissue. Observed molecular weight by Western-blot: 28-30 kDa

Preparation

This antibody was obtained by immunization of UBE2R2 recombinant protein (Accession Number: NM_017811). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116534(UBE2R2 Antibody) at dilution of 1:600