

# Human SCLY / Selenocysteine Lyase (His Tag) recombinant protein



Catalog Number: 500934

## General Information

### Protein Construction

A DNA sequence encoding the human SCLY (AAH07891.1) (Met1-Ala445) was expressed with a C-terminal polyhistidine tag.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Met

### Molecular Mass

The secreted recombinant human SCLY consists of 455 amino acids and predicts a molecular mass of 49.5 KDa. The apparent molecular mass of the protein is approximately 46 KDa in SDS-PAGE

under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human SCLY / Selenocysteine Lyase Protein (His Tag) SDS-PAGE