

GEFT antibody

Catalog Number: 110934

Product name

GEFT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GEFT Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Transfected HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of GEFT recombinant protein (Accession Number: NM_133483). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110934(GEFT antibody) at dilution of 1:300



IP Result of anti-GEFT (IP:Catalog No:110934, 4ug; Detection:Catalog No:110934 1:400) with HEK-293 cells lysate 2000ug.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110934(GEFT Antibody) at dilution of 1:1000