

IRSp53 antibody

Catalog Number: 111934

Product name

IRSp53 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

IRSp53 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human brain tissue. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, human brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 53 kDa


Validations



Preparation

This antibody was obtained by immunization of IRSp53 recombinant protein (Accession Number: NM_001385140). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111934(BAIAP2 Antibody) at dilution of 1:600



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:111934(BAIAP2 antibody) at dilution of 1:100 (under 10x lens)

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions

Recommended Dilution:



IRSp53 antibody



Catalog Number: 111934

IP Result of anti-BAIAP2 (IP:Catalog No:111934,

4ug; Detection:Catalog No:111934 1:500) with mouse brain tissue lysate 2640ug.