

# TEAD4 antibody

Catalog Number: 115934

## Product name

TEAD4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TEAD4 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MDA-MB-453s cells, BGC-823 cells, MCF-7 cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of TEAD4 recombinant protein (Accession Number: NM\_003213). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:115934 (TEAD4 Antibody) at dilution of 1:600



IP Result of anti-TEAD4 (IP:Catalog No:115934, 3ug; Detection:Catalog No:115934 1:300) with MCF-7 cells lysate 2000ug.