

# Cevostamab reference antibody



Catalog Number: 176934

## Product name

Cevostamab reference antibody

## Specificity

Human FcRH5 / IRTA2 / CD307e

## Antibody description

Anti- FcRH5 / IRTA2 / CD307e Antibody (Cevostamab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti- FcRH5 / IRTA2 / CD307e Antibody (cevoStamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti- FcRH5 / IRTA2 / CD307e Antibody (cevoStamab) is 91.59%, determined by



## SEC-HPLC.

Bioactivity: ELISA

Immobilized human IRTA2 His at 2 µg/mL can bind Anti- FcRH5 / IRTA2 / CD307e Antibody (cevoStamab),  $EC_{50}=0.007032$  µg/mL.