

# Anti-FTO Rabbit antibody



Catalog Number: 168934

## Product name

Anti-FTO Rabbit antibody

## Specificity

Human

## Antibody description

FTO Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human FTO

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/50

## Validations



Western blot detection of FTO in HeLa, Jurkat cell lysates using FTO Rabbit mAb(1:500 diluted). Predicted band size: 58kDa. Observed band size: 58kDa.