

Cavia porcellus CTLA4 (His tag) recombinant protein



Catalog Number: 518934

General Information

Gene Name Synonym

Cytotoxic T-lymphocyte protein 4; Cytotoxic T-lymphocyte-associated antigen 4; CTLA-4; CD152; CTLA4

Protein Construction

A DNA sequence encoding the Cavia porcellus CTLA4 Ala37-Asp161 is expressed with a 6His tag at the C-terminus.

Organism

Cavia porcellus

Expression Host

Yeast

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of PBS, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Cavia porcellus CTLA4 (His tag) recombinant protein