

Mouse CCL17/TARC recombinant protein



Catalog Number: 501234

General Information

Protein Construction

A DNA sequence encoding the mouse CCL17 (AAD56601.1) (Met 1-Pro 93) was expressed and purified.

Organism

Mouse

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 21

Molecular Mass

The secreted pro form of recombinant mouse CCL17 consists of 73 amino acids and has a

predicted molecular mass of 8.2 kDa as estimated in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CCL17 / TARC Protein SDS-PAGE