

# NR2B antibody

Catalog Number: 113234

## Product name

NR2B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NR2B Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 166kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000834). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113234(GRIN2B antibody) at dilution of 1:400



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:113234(GRIN2B Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)