

CRTC1 Monoclonal Antibody



Catalog Number: 166234

Product name

CRTC1 Monoclonal Antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Monoclonal antibody to CRTC1

Preparation

Antigen: Recombinant protein of human CRTC1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications


WB, ICC

Dilutions

WB 1:1000 - 1:4000

ICC 1:50 - 1:200

Validations



Western blot - CRTC1 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using CRTC1 antibody at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.