

PDP2 antibody

Catalog Number: 113634

Product name

PDP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PDP2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, mouse kidney tissue. Positive IHC detected in mouse brain tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of PDP2 recombinant protein (Accession Number: NM_001329933). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113634(PDP2 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:113634(PDP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

PDP2 antibody



Catalog Number: 113634

No:113634(PDP2 Antibody) at dilution of 1:50 (under 40x lens)