

RFX3 antibody

Catalog Number: 114634

Product name

RFX3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RFX3 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 70-79 kDa

Preparation

This antibody was obtained by immunization of RFX3 recombinant protein (Accession Number: NM_001282117). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114634(RFX3 antibody) at dilution of 1:600



IP Result of anti-RFX3 (IP:Catalog No:114634, 4ug; Detection:Catalog No:114634 1:500) with HeLa cells lysate 1400ug.